

Model 3552-3553

# Slings for children

Dino (green) & Butterfly (pink)



SUREHANDS®



- Your child is fully and safely supported
- Comfortable sitting or lying position
- Bottom remains free
- Easy to fit
- Quick transfers to and from bed, bath, shower, toilet or wheelchair
- Made-to-measure sizes available


Specialists in innovative and reliable lift & care systems

[www.surehands.com](http://www.surehands.com)

# Slings for children

## ● Dimensions (inch / mm)

	XS	S
A	2.2 / 56	2.28 / 58
B	2.6 / 66	2.95 / 75
C	0.53 / 13.5	0.55 / 14
D	0.3 / 7.5	0.31 / 8
E	1.1 / 28	1.26 / 32
F	2.13 / 54	2.48 / 63
G	1.48 / 37.5	1.87 / 47.5
H	2.09 / 53	2.32 / 59
I	1.65 / 42	1.85 / 47



## ● Technical specifications

Materials ..... Polyester  
 Safe Working Load ..... 496 lbs / 225 kg

## ● Art. nrs.

3552 : Butterfly Pink  
 3553 : Dino Green

Local Representative

Lift & Care® is a registered trademark of  
 T.F. Herceg, Inc. dba SUREHANDS LIFT & CARE SYSTEMS.

© 04-04-2024 • Handi-Move • v1.2